



GenInst Studio Configuration Guide



Choose a software subscription to generate custom instruments and a hardware license to deploy them on your Moku.

GenInst Studio software subscription tiers

Select the GenInst Studio tier that meets your needs.

Base \$0*

Includes:

GenInst Studio Access

Build instruments for any Moku device

Up to 10 instruments/year

Build 1 instrument at a time

Limited execution, synthesis, routing times

Data may be used for product improvement

*Requires Moku registration

Premium \$5,000 USD

Everything in Base, plus:

Up to 100 instrument builds per year

Build 2 instruments in parallel

Download generated HDL code and testbench

Export projects to Moku Cloud Compile

Extended execution time

Interject updates during generation

Opt out of data usage for product improvement

Paid per user, per year

Don't have a Moku, but want to try GenInst Studio? [Access a 14-day trial](#) to build up to 3 instruments. Other Base tier limitations apply.

Moku device hardware license

To run your custom capability on a Moku device, select GenInst from the instrument menu. GenInst is available on all Moku devices free of charge.

GenInst

Deploy and run GenInst Studio builds for any Moku device.

| \$0, included in Base Bundle

Maximum Sample Rates

Moku:Delta			Moku:Pro	Moku:Lab		Moku:Go	
3-slot	5-slot	8-slot	4-slot	2-slot	3-slot	2-Slot	3-slot
312.5 MSa/s			312.5 MSa/s	125 MSa/s		31.25 MSa/s	

If you have any questions or need help configuring, please reach out to support@liquidinstruments.com.